



ES256H



ES256H GENERAL DISINFECTANT CONCENTRATE

- Phosphate-free formulation
- Neutral pH will not dull finishes with repeated use

	ES256H
APPEARANCE	Clear blue liquid
ODOR	Lemon
SOLUBILITY	Soluble
RELATIVE DENSITY	~1.011
pH- VALUE	6.8-7.8
SHELF LIFE	3 years

SIZES	
ES256-CS	4.73 L (1.25 US gallon)
ES256-4CS2	1.89 L (0.50 U.S. gallon)

BENEFITS

- One step neutral disinfectant that is effective against a broad spectrum of bacteria, is virucidal* (including HIV-1, HIV-2, HBV and HCV) and inhibits the growth of mold and mildew and their odors when used as directed
- Cleans, shines, deodorizes and disinfects all hard non-porous surfaces listed on the label
- Neutral pH multi surface cleaner, deodorizer and disinfectant

ACTIVE INGREDIENTS GUARANTEE

Octyl decyl dimethyl ammonium chloride	6.510%
Diocetyl dimethyl ammonium chloride	3.255%
Didecyl dimethyl ammonium chloride	3.255%
Alkyl (50% C14, 40% C12, 10% C16) dimethyl benzyl ammonium chloride	8.680%

USES:

Use on floors, walls, toilet bowls, urinals and other hard non-porous surfaces such as sinks, tables and counter tops made of Stainless Steel, Metal, Glazed Porcelain. Other areas of use include:

- Farms
- Florists
- Hospitals/Nursing Homes
- Hotels
- Kitchens
- Medical, dental and veterinary offices
- Schools

Can be used with mop and bucket, trigger sprayer, sponge soaking or by auto scrubber.

DILUTIONS

USE	RATIO
Disinfectant/Cleaner	1:256

CAUTION: Refer to label for complete instructions, dilutions, safety requirements and kill claims

DISPOSAL

Dispose of used/unwanted product according to state and local regulations. Offer container for recycling.

PERSONAL PROTECTION

Please refer to product label and SDS.

DISINFECTANTS