



## ES65 HYDROGEN PEROXIDE DISINFECTANT CLEANER

- Oxidizing multi-purpose degreaser
- Fragrance and Dye Free
- 5 minute registered dwell time
- Effective against Herpes Simplex Type 2, Influenza Type A2 and HIV-1

	ES65
APPEARANCE	Clear liquid
ODOUR	None
SOLUBILITY	Soluble
RELATIVE DENSITY	1.01
pH-VALUE	5.4 - 5.8
VISCOSITY	Water thin
SHELF LIFE	1 year

SIZES	
ES65-CS	4.73 L (1.25 U.S. gallon)
ES65-1S	2 L (0.52 U.S. gallon)
ES65-4	3.78 L (1.0 U.S. gallon)

### BENEFITS

- Oxidizing multi-purpose degreaser
- Effective against Herpes Simplex Type 2, Influenza Type A2 & HIV-1
- 5 minute dwell time
- One product program with multiple dilutions
- Safer, lockable cabinetry

### ACTIVE INGREDIENT

Hydrogen Peroxide.....3.95%  
Benzalkonium Chloride .....2.50%

### USES

- Hospitals
- Medical and dental
- Clinics
- Nursing homes
- Schools
- Offices
- Food Service Establishments
- Locker rooms
- Restrooms
- Playrooms
- Daycares
- Veterinary clinics
- Hotels
- Exercise facilities

### DILUTIONS

USE	RATIO	PER GALLON OF WATER	PER LITRE OF WATER
Disinfectant	1:12	10.7 oz	317 mL
Cleaner	1:127	1 oz	7.8 mL

**NOTE:** Improper use of dilutions may cause this product not to work properly and may increase the risk of health effects.

### DIRECTIONS

- Use with cold water
- Remove gross filth and heavy soil deposits
- A pre-cleaning step is required and recommended for heavily soiled areas
- Apply solution with a spray bottle, cloth, sponge or mop to wet thoroughly
- Allow to remain wet for 5 minutes and then let air dry
- Prepare a fresh solution for each use

### DISPOSAL

Dispose of used/unwanted product according to state and local regulations. Offer container for recycling.

### PERSONAL PROTECTION

Please refer to product label and SDS.